

Drug Delivery Using Nanoparticles

Nicomp® DLS

DRUG DELIVERY USING NANOPARTICLES

Considerable research and development is devoted to improving drug delivery through the use of nanoparticles. Although the definition of nanoparticle has become murky, typically they are defined being in the size range of 100 nm and below, Dynamic light scattering (DLS) is the preferred method in this size range and the Nicomp® is often used in this field of research. This application note summarizes how the Nicomp can be used to study and control nanoparticles used for drug delivery.

INTRODUCTION

Both ISO/TS 27687¹ and ASTM E2456² define nanoparticles as being in the size range of 100 nm and below, making this the most widely used classification. Less strict interpretations have extended the upper size range for both scientific, and other reasons. Now many nanomaterials greater than 100 nm in size are commonly called nanoparticles. The motivations for developing drug products in this size range include improved dissolution/bioavailability, targeting, circulation time in the system, and area under the curve (AUC) pharmacokinetics.

Many of these drug products are developed to enhance targeting. A passive targeting approach increases the circulation time by reducing the size and cloaking the nanoparticle with a coating such as polyethylene glycol (PEG). An active targeting approach modifies the surface of the nanoparticle to seek and adhere to specific parts of the body, such as cancer tumors, while avoiding healthy tissue. Cell specific ligands on the surface of the nanoparticle can be added to bind specifically to complementary receptors.

The Nicomp DLS system (Figure 1) is an ideal instrument to measure both the size and zeta potential (surface charge) of nanoparticles used for drug delivery.



Figure 1. Nicomp DLS instrument

TYPES OF NANOPARTICLES

Nanocrystals

Active pharmaceutical ingredients (APIs) are often crystalline. Hydrophobic crystals can be difficult to formulate to be delivered in a hydrophilic carrier mechanism. By reducing the size to the nanocrystal range, a nanosuspension can improve the bioavailability of drugs, where the dissolution velocity is the rate limiting step, such as poorly water soluble drugs.³ These nanocrystal often need to be stabilized using surfactants or polymers including during processing. A decrease in particle size increases dissolution rate by both increasing the surface area A (Figure 2) and the saturation solubility C_s .




Figure 2. Surface enlargement with size reduction

The Noyes-Whitney equation (Equation 1) shows how an increase in both A and C_s will affect the dissolution rate dC/dt .

$$\frac{dC}{dt} = \frac{DA}{Vh} (C_s - C_x) \dots \dots \dots \text{(Equation 1)}$$

Where

- dC/dt = dissolution rate
- D = diffusion coefficient
- A = surface area
- C_s = concentration at boundary layer
- C_x = concentration API @ given time
- V = volume dissolution medium
- h = height of boundary layer

Lipid-based liquid crystalline nanoparticles (LCNPs) are another delivery system capable of increasing bioavailability of both hydrophobic and hydrophilic drugs. These are self-assembled structures prepared by high shear energy dispersing of a nonlamellar liquid crystalline matrix into the water phase. The particle size of the LCNPs is an important physiochemical property requiring proper analysis and control. The Nicomp DLS system has been successfully used to determine both the mean size, and presence of aggregates in LCNP dispersions.⁴ Paclitaxel was loaded into a LCNP dispersion and analyzed by the Nicomp DLS system and TEM, see Figure 3.




Figure 3. Nicomp and TEM results for an LCNP dispersion, copied with rights from⁴

The TEM image indicates a bimodal particle size distribution of smaller near, 25 nm particles, plus larger particles on the scale of 100 nm. The upper Nicomp result is the Gaussian intensity distribution mean forcing the entire distribution into one peak. The lower Nicomp result utilizes the proprietary Nicomp non-negative least squares algorithm to

report a higher resolution and more accurate description of the bimodal nature of the actual particle size distribution. This highlights a main advantage of the Nicomp DLS system – the ability to resolve multi-modal distributions even at concentrations as low as 0.2 mg/mL.⁵

Micelles

Another potential drug delivery system for increasing the solubilization of hydrophobic drugs is polymeric micelles.⁶ Micelles are formed when the concentration of the polymer, in solution, exceeds a certain threshold concentration known as the critical micellar concentration (CMC). Polymeric micelles are core-shell nanostructures synthesized from amphiphilic block copolymers. Micelles have the advantages of being very small in size (10 – 100 nm), improving passive targeting to solid tumors. By modifying the surface with ligands polymeric micelles can become capable of site-specific drug delivery.

The Nicomp DLS system has been used for particle size measurements in many micelle based research projects.⁷⁻¹¹ In one study,¹² polymeric micelles were formed using copolymers polycaprolactone (PCL) and polyethylene glycol (PEG). Docetaxel (DTX) was used as the model drug and the surface was modified with a small molecular ligand of prostate specific membrane antigen (PSMA). Figure 4 shows the self-assembly of the micelles and the endocytosis process of the drug loaded final structure.




Figure 4. Preparation and endocytosis of DTX loaded polymeric micelles targeted to PSMA¹²

The particle size by the Nicomp DLS system, and TEM of two samples used in this study, is shown in Figure 5. The data for non-targeted micelles are shown on the left, and the targeted on the right. The DLS data appears slightly larger than the TEM images, possibly due to shrinkage of the PEG shell induced by water evaporation before TEM analysis.




Figure 5. Size of non-targeted (upper) and targeted (lower) polymeric micelles by DLS and TEM¹²

Liposomes

Liposomes are bilayer vesicles routinely used in the pharmaceutical industry as a drug delivery system for transport of chemotherapeutic drugs to the tumor area. They are composed of phospholipids that have a polar end attached to a nonpolar chain that self-assemble into bilayer vesicles with the polar ends facing the aqueous medium and nonpolar ends forming a bilayer. In pharmaceutical applications, the active pharmaceutical ingredient (API) is usually incorporated into the liposome either into the hydrophilic pocket or sandwiched between the bilayers depending on the hydrophilicity of the API, see Figure 6. Surface modification is common for targeted delivery.




Figure 6. Complex liposome structure

Monitoring the particle size while processing liposomes is critical and the Nicomp DLS system is frequently used for this application.¹³⁻²⁰ In one internal Entegris study, liposomes were created using a formulation of 3:1:1 HSPC, cholesterol and mPEG-DSPE. The sample was first mixed by rotor stator at 7200 rpm for 10 minutes, then passed through a Microfluidizer²¹ at 25,000 psi using a Y chamber to create the liposomes. The samples were processed 1, 3, 5, and 10 passes through the microfluidizer. An image of the pre-mix and processed samples (left to right) is shown in Figure 7.




Figure 7. Pre-mix, 1, 3, 5, and 10 passes

The liposome samples were analyzed on both the Nicomp DLS system and the AccuSizer® single particle optical sizing (SPOS) system. DLS was used to determine the reduction of the intensity mean size during processing, while the AccuSizer (LE sensor range 0.5 – 400 μm) was used to quantify the presence of larger particle tails of the distribution. The Nicomp DLS results are shown in Figure 8, and the AccuSizer SPOS results are shown in Figure 9.




Figure 8. Nicomp DLS results from right to left; premix, 1, 3, 5, and 10




Figure 9. AccuSizer SPOS results right to left; premix, 1, 3, 5, and 10 passes

Using both DLS to determine mean size and SPOS to quantify the presence and concentration of tails is common in many industries and is an integral part of USP <729> Globule-size distribution in lipid injectable emulsions.²²

Online DLS for process monitoring

While the vast majority of DLS measurements are made in the laboratory, Entegris has installed several systems in customer manufacturing operations that track particle size during production runs.²³ These systems

have been used to monitor high-pressure homogenization processes used during the manufacture of nanoparticles for drug delivery. The at-line system removes a sample from the process, dilutes the sample to avoid multiple scattering effects, measures the sample, and then repeats the procedure (see Figure 10). The complete measurement cycle is approximately two minutes, providing continuous particle size information to the process engineers monitoring the manufacturing operation.




Figure 10. Online DLS system schematic

Figure 11 shows online DLS results as a function of pressure downstream of a high pressure homogenizer. The goal was to determine the optimum pressure to keep the particle size very close to 100 nm in size. After the optimum pressure (~10 k psi) was determined the online DLS system was used to assure the complete batch was manufactured within specification.




Figure 11. Pressure vs. particle size in process DLS results

CONCLUSIONS

The Nicomp DLS system is widely used for particle size and zeta potential analysis of drug delivery systems in the nano scale in research,²⁴⁻³⁹ quality release testing and in process monitoring. The AccuSizer SPOS provides a complementary technique for determining the concentration of larger particles that could indicate instability or non-optimized formulation or process conditions.

References

- ¹ ISO/TS 27687, Nanotechnologies — Terminology and definitions for nanoobjects — Nanoparticle, nanofibre and nanoplate, available at <https://www.iso.org/obp/ui/#iso:std:iso:ts:27687:ed-1:v2:en>
- ² ASTM E2456, Standard Terminology Relating to Nanotechnology, available at Standards/E2456.htm
- ³ Jens-Uwe et al., Nanocrystal technology, drug delivery and clinical applications, International Journal of Nanomedicine 2008:3(3) 295–309
- ⁴ Zeng et al., Lipid-based liquid crystalline nanoparticles as oral drug delivery vehicles for poorly water-soluble drugs International Journal of Nanomedicine 2012:7
- ⁵ Scomparin et al., Novel folated and non-folated pullulan bioconjugates for anticancer drug delivery European Journal of Pharmaceutical Sciences 42 (2011) 547–558
- ⁶ Cory et al., Polymeric Micelles for Drug Delivery, Curr Pharm Des. 2006;12(36):4669-84
- ⁷ Koizumi et al., Novel SN 38 Incorporating Polymeric Micelles, NK012 Eradicate Vascular Endothelial Growth Factor Secreting Bulky Tumors, Cancer Res 2006; 66: (20) with Nicomp data
- ⁸ Song et al., Self-assembled micelles of novel amphiphilic copolymer cholesterol-coupled F68 containing cabazitaxel as a drug delivery system, Int J Nanomedicine. 2014; 9: 2307–2317.
- ⁹ Wang, Pharmacokinetics and Biodistribution of Paclitaxel-loaded Pluronic P105/L101 Mixed Polymeric Micelles, Pharmaceutical Society of Japan, 128(6), 2008
- ¹⁰ Bachar et al., Development and characterization of a novel drug nanocarrier for oral delivery, based on self-assembled β -casein micelles, Journal of Controlled Release, Volume 160, Issue 2, 10 June 2012
- ¹¹ Jiang et al., Poly(aspartic acid) derivatives as polymeric micelle drug delivery systems J Appl Polym Sci 101: 2871–2878, 2006
- ¹² Jin et al., PSMA Ligand Conjugated PCL-PEG Polymeric Micelles Targeted to Prostate Cancer Cells, PLoS ONE 9(11): e112200. doi:10.1371/journal.pone.0112200
- ¹³ Zidan et al., Near-Infrared Investigations of Novel Anti-HIV Tenofovir Liposomes, The AAPS Journal, Vol. 12, No. 2, June 2010
- ¹⁴ Wong et al., A New Polymer-Lipid Hybrid Nanoparticle System Increases Cytotoxicity of Doxorubicin Against Multidrug-Resistant Human Breast Cancer Cells, Pharmaceutical Research, Vol. 23, No. 7, July 2006
- ¹⁵ Zhang et al., The cargo of CRPPR-conjugated liposomes crosses the intact murine cardiac endothelium, J Control Release, 2012 October 10; 163(1)
- ¹⁶ Guan et al., Enhanced oral bioavailability of cyclosporine A by liposomes containing a bile salt, International Journal of Nanomedicine 2011:6
- ¹⁷ Ando et al., Reactivity of IgM antibodies elicited by PEGylated liposomes or PEGylated lipoplexes against auto and foreign antigens, Journal of Controlled Release, Volume 270, 28 January 2018
- ¹⁸ Johnston et al., Characterization of the drug retention and pharmacokinetic properties of liposomal nanoparticles containing dihydrosphingomyelin, Biochimica et Biophysica Acta 1768 (2007)
- ¹⁹ Cipolla et al., Modifying the Release Properties of Liposomes Toward Personalized Medicine, Journal of Pharmaceutical Sciences 103:1851–1862, 2014
- ²⁰ El-Ridy et al., Liposomal Encapsulation of Amikacin Sulphate for Optimizing Its Efficacy and Safety, BJPR, 5(2): 98-116, 2015
- ²¹ Entegris Application Note Size Reduction by a Microfluidizer,
- ²² Entegris Application Note USP 729 Testing
- ²³ Entegris Application Note Nanoparticles for Drug Delivery
- ²⁴ Wong et al., A New Polymer-Lipid Hybrid Nanoparticle System Increases Cytotoxicity of Doxorubicin Against Multidrug-Resistant Human Breast Cancer Cells, Pharmaceutical Research, Vol. 23, No. 7, July 2006
- ²⁵ Martins et al., Brain delivery of camptothecin by means of solid lipid nanoparticles: Formulation design, in vitro and in vivo studies, International Journal of Pharmaceutics 439 (2012) 49– 62
- ²⁶ Podaralla et al., Influence of Formulation Factors on the Preparation of Zein Nanoparticles, AAPS PharmSciTech, Vol. 13, No. 3, September 2012
- ²⁷ Chertok et al., Iron oxide nanoparticles as a drug delivery vehicle for MRI monitored magnetic targeting of brain tumors, Biomaterials, Volume 29, Issue 4, February 2008
- ²⁸ Songa et al., Formulation and characterization of biodegradable nanoparticles for intravascular local drug delivery, Journal of Controlled Release, Volume 43, Issues 2–3, 18 January 1997
- ²⁹ Jain et al., Magnetic nanoparticles with dual functional properties: Drug delivery and magnetic resonance imaging, Biomaterials, Volume 29, Issue 29, October 2008
- ³⁰ Guo et al., Aptamer-functionalized PEG–PLGA nanoparticles for enhanced anti-glioma drug delivery, Biomaterials Volume 32, Issue 31, November 2011
- ³¹ Nguone et al., Accumulating nanoparticles by EPR: A route of no return, Journal of Controlled Release Volume 238, 28 September 2016Menzel et al., In vivo evaluation of an oral self-emulsifying drug delivery system (SEDDS) for exenatide, Journal of Controlled Release, Volume 277, 10 May 2018
- ³² Dorati et al., Gentamicin Sulfate PEG-PLGA/PLGA-H Nanoparticles: Screening Design and Antimicrobial Effect Evaluation toward Clinic Bacterial Isolates, Nanomaterials 2018, 8, 37
- ³³ Xu et al., The performance of docetaxel-loaded solid lipid

nanoparticles targeted to hepatocellular carcinoma, *Biomaterials* 30 (2009) 226–232

³⁴ Piao et al., *Human serum albumin-coated lipid nano-particles for delivery of siRNA to breast cancer*, *Nanotechnology, Biology, and Medicine* 9 (2013)

³⁵ Andersen et al., *Chitosan-Based Nanomedicine to Fight Genital Candida Infections: Chitosomes*, *Mar. Drugs* 2017, 15, 64

³⁶ Kou et al., *Preparation and characterization of the Adriamycinloaded amphiphilic chitosan nanoparticles and their application in the treatment of liver cancer*, *Oncology Letters* 17: 7833-7841, 2017

³⁷ Kuang et al., *Dual Functional Peptide-Driven Nano-particles for Highly Efficient Glioma-Targeting and Drug Codelivery*, *Molecular Pharmaceutics*, April, 2016

³⁸ Cooper et al., *Formulation and in vitro evaluation of niacin-loaded nanoparticles to reduce prostaglandin mediated vasodilatory flushing*, *European Review for Medical and Pharmacological Sciences*, 2015; 19: 3977-3988

³⁹ Martins et al., *Physicochemical properties of lipid nanoparticles: Effect of lipid and surfactant composition*, *Drug Development and Industrial Pharmacy* 37(7):815-24

FOR MORE INFORMATION

Please call your Regional Customer Service Center today to learn what Entegris can do for you. Visit entegris.com and select the [Contact Us](#) link to find the customer service center nearest you.

TERMS AND CONDITIONS OF SALE

All purchases are subject to Entegris' Terms and Conditions of Sale. To view and print this information, visit entegris.com and select the [Terms & Conditions](#) link in the footer.



Corporate Headquarters

129 Concord Road
Billerica, MA 01821
USA

Customer Service

Tel +1 952 556 4181
Fax +1 952 556 8022
Toll Free 800 394 4083

Entegris®, the Entegris Rings Design®, and other product names are trademarks of Entegris, Inc. as listed on entegris.com/trademarks. All third-party product names, logos, and company names are trademarks or registered trademarks of their respective owners. Use of them does not imply any affiliation, sponsorship, or endorsement by the trademark owner.

©2018-2019 Entegris, Inc. | All rights reserved. | Printed in the USA | 7130-10549TAN-0819