

Development of Alternative BCDP Reagents for the Synthesis of 5'-Terminal Phosphate Oligonucleotides: Application in the Synthesis of SmartCap™, a Novel mRNA Capping Reagent

Smita Patnaik, Scott Laneman, Jonathan Dube – Entegris

Dong Hun Jung, Jaemun Jung, Mun Sik Seo, Jun Young Choi – ST PHARM


ABSTRACT

BCDP (Bis(2-cyanoethyl) diisopropylphosphoramidite) is a phosphitylating agent that facilitates the formation of phosphite triester bonds between nucleosides. These bonds can be subsequently oxidized or sulfurized to produce a terminal phosphate or thiophosphate at the 3' or 5' hydroxyl group of various oligonucleotides. Terminal phosphate or thiophosphate oligonucleotides serve as versatile building blocks for chemical or enzymatic conjugation reactions, enabling the synthesis of longer nucleotide sequences.

In this technical note, we will address the safety concerns associated with the use of BCDP reagent and explore the development of alternative phosphitylating reagents. Additionally, we will discuss its application in the synthesis of SmartCap™, a novel mRNA capping reagent.

INTRODUCTION


A General Approach to Terminal Modification via BCDP Phosphitylation



BCDP is a phosphoramidite commonly used in oligonucleotide synthesis, but it requires careful handling due to concerns related to stability, exothermicity, and toxicity.

Risk information about Bis(2-cyanoethyl)-N,N-diisopropylphosphoramidite. 'BCDP'

Exothermicity	- UN class 1 explosive - Flammable liquid and vapor
Stability	- Recommended storage temperature -20°C - Handle and store under inert gas - At room temperature, the purity decreased from 96.88% to 8.22% within 30 days
Toxicity	- Class 6.1 in DOT - May cause an allergic skin reaction



Alternative P-reagents were introduced at this stage of the process.

Thermal Stability (ARC and DSC) of P-Reagent Substitutes


ARC Testing Comparison of Phosphoramidites

Structure				
Compound	BCDP ~98% Purity	Substitute-1 96% purity	Substitute-2 ~99% purity	Substitute-3
Earliest onset temperature	75.9°C	170.9°C	170.5°C	130.7°C
Earliest onset energy	-473.0 J/g	-490.5 J/g	-142.9 J/g	-224.4 J/g
Number of additional onsets	3	2	2	2


Table continues on next page

Highest energy onset	-772.3 J/g (190°C)	-409.5 J/g (170.9°C)	>-158.5 J/g (330.1°C)	-224.4 J/g (130.7°C)
Time to maximum rate (TMR) estimates	0.9 hr at 70°C 24 hr at 48°C 70 hr at 41°C	8 hr at 156°C 24 hr at 143°C 48 hr at 128°C	0.5 hr at 176°C 24 hr at 170°C 48 hr at 169°C	0.8 hr at 130°C 24 hr at 118.5°C 71 hr at 115°C
Total heat released	-1,271.3 J/g	-616.6 J/g	-301.4 J/g	-348.4 J/g

Temperature and Pressure vs. Time BCDP
Accelerating Rate Calorimetry



Temperature and Pressure vs. Time of Substitute-3
Accelerating Rate Calorimetry



DSC Testing Comparison of Phosphoramidites


SAMPLE	EVENT	START TEMP.	PEAK TEMP.	END TEMP.	ENERGY
BCDP	Exotherm	66.7°C	129.4°C	194.0°C	-542.7 J/g
	Exotherm	194.1°C	284.6°C	382.6°C	-165.4 J/g
Substitute-3	Endotherm	35.3°C	47.2°C	55.6°C	97.74 J/g
	Exotherm	190.6°C	206.6°C	217.5°C	-232.45 J/g
	Endotherm	323.3°C	362.6°C	390.0°C	317.53 J/g

BCDP is a candidate for classification as a UN Class 1 explosive substance.

However, based on the thermograms for 'substitute-3', this chemical is not a candidate for classification as a UN Class 1 explosive substance or a UN Class 4, Division 4.1 self-reactive substance.


APPLICATION OF ALTERNATIVE P-REAGENTS IN THE SYNTHESIS OF SMARTCAP

Step 1 – Comparative Phosphitylation Results Using Various Phosphoramidites




PHOSPHORAMIDITES	PHOSPHINE ACTIVATOR	IPC RESULT (HPLC AREA %)	
		RESIDUAL B-3	RXN TIME (HR)
BCDP	1H-Tetrazole	1.40 ~ 4.50%	1
Substitute-1	1H-Tetrazole	0.42%	1
Substitute-2	1H-Tetrazole	~ 0.50%	~ 2
Substitute-3	1H-Tetrazole	~ 0.04%	~ 3
	5-(Benzylthio)-1H-tetrazole	0.06%	1
	5-(Ethylthio)-1H-tetrazole	0.17%	1
	4,5-dicyanoimidazole	0.45%	1

Step 2 – Comparative Oxidation Results Using Various Phosphoramidites



PHOSPHORAMIDITES	OXIDATION REAGENT	IPC RESULT (HPLC AREA %)		ISOLATION YIELD
		RESIDUAL B-2	RXN TIME (HR)	
BCDP	I ₂ , (THF : PW : Py = 66 : 33 : 1)	~ 0.50%	~ 1	63 ~ 72%
Substitute-1	H ₂ O ₂	0.37%	~ 1	68%
Substitute-2	H ₂ O ₂	0.35%	~ 1	30%
Substitute-3	I ₂ , (THF : PW : Py = 66 : 33 : 1)	N.D	~ 1	67%

Step 3 – Comparative Deprotection Results Using Various Phosphoramidites



PHOSPHORAMIDITES	IPC RESULT (HPLC AREA %)		ISOLATION YIELD
	RESIDUAL B-1	RXN TIME (HR)	
BCDP	N.D ~ 1.3%	24	43 ~ 53%
Substitute-3	1.4%	48	67%

Using the same process, deprotection with Substitute-3 enabled the synthesis of the dinucleotide with a higher yield compared to BCDP.

CONCLUSION

In the manufacturing process of dinucleotides for SmartCap synthesis, BCDP is used as the phosphitylating reagent. However, BCDP is susceptible to oxidation to P(V) during storage and inherently exhibits explosive properties.

In contrast, ARC and DSC analyses confirmed that the three alternative reagents exhibit reduced hazardous characteristics. Furthermore, it was experimentally demonstrated that all three alternatives can effectively replace BCDP in the synthesis of dinucleotides required for SmartCap production.

ACKNOWLEDGEMENTS

This work was supported by the Technology Innovation Program (Grant #: 20018209) funded by the Ministry of Trade, Industry & Energy (MOTIE, South Korea).

The alternative reagent for the BCDP replacement study was provided by Entegris.

FOR MORE INFORMATION

Please call your Regional Customer Service Center today to learn what Entegris can do for you. Visit [entegris.com](https://www.entegris.com) and select the [Contact Us](#) link to find the customer service center nearest you.

TERMS AND CONDITIONS OF SALE

All purchases are subject to Entegris' Terms and Conditions of Sale. To view and print this information, visit [entegris.com](https://www.entegris.com) and select the [Terms & Conditions](#) link in the footer.



Corporate Headquarters

129 Concord Road
Billerica, MA 01821
USA

Customer Service

Tel +1 952 556 4181
Fax +1 952 556 8022
Toll Free 800 394 4083

Entegris®, the Entegris Rings Design®, and other product names are trademarks of Entegris, Inc. as listed on [entegris.com/trademarks](https://www.entegris.com/trademarks). All third-party product names, logos, and company names are trademarks or registered trademarks of their respective owners. Use of them does not imply any affiliation, sponsorship, or endorsement by the trademark owner.

©2026 Entegris, Inc. | All rights reserved. | Printed in the USA | 9000-14681ENT-0726