

# Nicomp® Online Nanoparticle Size Analyzer

Now measure nanoparticles 1 nm – 1 μm continuously with the Nicomp Online system

Dynamic light scattering (DLS) is the most prevalent size analysis technique for nanoparticles. Although the vast majority of DLS size measurements are performed in the laboratory, Entegris has been supplying online DLS system solutions for over 25 years.<sup>1,2</sup> The nanotechnology revolution has created higher demand for online DLS systems and the Nicomp® Online is a proven industrial solution.

**Operation:** The Nicomp Online schematic is shown in Figure 1. The pressurized sample line enters a variable diluter where the diluent automatically adjusts the concentration to the optimum range for the measurement. The diluted sample enters the scattering cell where the DLS measurement is performed. Once the measurement is completed, the sample is pumped to drain, the system is cleaned and another measurement is initiated.

**Process integration is easy:** The optics units measure 10"x12"x6" as shown in Figure 2. The sample and diluent enter the system via 0.25" tubing. The inlet sample pressure can be 10 – 1500 psi. An electronics control box (not shown) connects to the optics/fluidics enclosures seen in Figure 2 and can be located 1 – 5 m from the optics/fluidics. A standard process environmentally hardened computer controls the Nicomp Online system.

**Results:** Figure 3. Shows data collected at a customer site<sup>3</sup> with mean particle size on the left y-axis and homogenizer pressure on the right y-axis. Results were generated ~ every 3 minutes, provided the ability to optimize operating conditions to achieve the desired size result.




Figure 1. Optics/fluidics schematic



Figure 2. Optics/fluidics module




Figure 3. Size vs. pressure results

FOR INTERNAL PRINT APPROVAL ONLY • FOR INTERNAL PRINT APPROVAL ONLY

## SPECIFICATIONS\*

<b>Principle</b>	Dynamic light scattering (DLS), 90° data collection angle
<b>Size range</b>	1 nm – 1 µm, depending on sample density
<b>Analysis algorithms</b>	Gaussian and Nicomp high resolution multi-modal with ability to resolve peaks separated by 1.5X in size (220 and 340 nm)
<b>Detectors</b>	PMT standard, optional 7X gain APD
<b>Light adjustment</b>	Light intensity automatically controlled using a neutral density filter
<b>Temperature range</b>	10° – 50°C
<b>Physical</b>	Electronics enclosure: 10"w x 12"h x 6"d Fluidics/optics enclosure: 10"w x 12"h x 6"d
<b>Power</b>	100 – 120 VAC, 60 Hz or 220 – 240 VAC, 50 Hz

\* The Nicomp Online system includes three components: the optics/fluidics enclosure, the electronics enclosure and a computer. The specifications presented in this flyer are typical, but custom systems built for specific processes and requirements are possible.

## References

- <sup>1</sup> Kourti, T., et. al., Online Particle Size Determination during Latex Production using Dynamic Light Scattering, American Chemical Society, 1991
- <sup>2</sup> Nicoli, D., et. al., Online submicron particle sizing : Carrying the analytical laboratory to the process facility, American Laboratory, October, 1990
- <sup>3</sup> Data courtesy of Bind Therapeutics, Cambridge, MA

## FOR MORE INFORMATION

Please call your Regional Customer Service Center today to learn what Entegris can do for you. Visit [entegris.com](http://entegris.com) and select the [Contact Us](#) link to find the customer service center nearest you.

## TERMS AND CONDITIONS OF SALE

All purchases are subject to Entegris' Terms and Conditions of Sale. To view and print this information, visit [entegris.com](http://entegris.com) and select the [Terms & Conditions](#) link in the footer.



### Corporate Headquarters

129 Concord Road  
Billerica, MA 01821  
USA

### Customer Service

Tel +1 952 556 4181  
Fax +1 952 556 8022  
Toll Free 800 394 4083

Entegris®, the Entegris Rings Design®, and other product names are trademarks of Entegris, Inc. as listed on [entegris.com/trademarks](http://entegris.com/trademarks). All third-party product names, logos, and company names are trademarks or registered trademarks of their respective owners. Use of them does not imply any affiliation, sponsorship, or endorsement by the trademark owner.

©2018-2019 Entegris, Inc. | All rights reserved. | Printed in the USA | 7130-10471TAN-0519