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### Gaussian Summary

|   |                              |
|---|------------------------------|
| Mean diameter                                   | 132.5 nm                     |
| Standard deviation                              | 44.0 nm (33.2%)              |
| Normalized standard deviation (Coeff. of Var'n) | 0.332                        |
| Variance (P.I.)                                 | 0.110                        |
| Chi squared                                     | 0.577                        |
| Baseline Adj.                                   | 0.000%                       |
| Z-Avg. Diff. Coeff                              | 3.51E-008 cm <sup>2</sup> /s |

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### Gaussian Summary

|   |                              |
|---|------------------------------|
| Mean diameter                                   | 122.8 nm                     |
| Standard deviation                              | 16.2 nm (13.2%)              |
| Normalized standard deviation (Coeff. of Var'n) | 0.132                        |
| Variance (P.I.)                                 | 0.017                        |
| Chi squared                                     | 0.173                        |
| Baseline Adj.                                   | 0.078%                       |
| Z-Avg. Diff. Coeff                              | 3.78E-008 cm <sup>2</sup> /s |

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### Gaussian Summary

|   |                              |
|---|------------------------------|
| Mean diameter                                   | 139.3 nm                     |
| Standard deviation                              | 30.8 nm (22.1%)              |
| Normalized standard deviation (Coeff. of Var'n) | 0.221                        |
| Variance (P.I.)                                 | 0.049                        |
| Chi squared                                     | 2.059                        |
| Baseline Adj.                                   | 0.000%                       |
| Z-Avg. Diff. Coeff                              | 3.34E-008 cm <sup>2</sup> /s |

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- <sup>1</sup> Aoki, N. & Hashimoto, M., Hashimoto Electronic Industry CO. and Yoshimura, T. Liposome Engineering Laboratory, *Measurement of Liposome Size Distribution Using Nicomp 380 and AccuSizer 780 AD*, presentation, July 2013
- <sup>2</sup> Thanks to Elizabeth Ingham, Azadeh Kheiriloom and Jai Seo from the Dr. Katherine Ferrara Lab in the Department of Biomedical Engineering at UC Davis for sharing these data and helping to create this document.
- <sup>3</sup> Mahakian, L. et.al., *Comparison of PET Imaging with <sup>64</sup>Cu-Liposomes and <sup>18</sup>F-FDG in the 7,12-Dimethylbenz[a] anthracene (DMBA)-Induced Hamster Buccal Pouch Model of Oral Dysplasia and Squamous Cell Carcinoma*, Mol Imaging Biol (2013)
- <sup>4</sup> Kheiriloom, A. et.al., *Complete regression of local cancer using temperature-sensitive liposomes combined with ultrasound-mediated hyperthermia*, Journal of Controlled Release, 172 (2013)
- <sup>5</sup> Seo, J. et.al., *The pharmacokinetics of Zr-89 labeled liposomes over extended periods in a murine tumor model*, Nuclear Medicine and Biology 42 (2015)

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©2018-2021 Entegris, Inc. | All rights reserved. | 7127-10535ENT-0421JPN